



**Alpen**<sup>®</sup>  
**POLISHER INSTRUMENTS**

## Provides a smoother polishing experience

- An ideal offering of Silicone and Diamond Impregnated Polishers
- Provides the clinician with a choice of 1, 2 or 3-step process
- Polisher designed to handle the unique characteristics of porcelain veneers, inlays and onlays – chairside and in the lab
- Allows for fine contouring to a mirror shine
- Long lasting shape integrity

**coltène**   
**whaledent**<sup>®</sup>




## Polishers that work with delicate materials, producing astounding results

- Alpen Polishers produce a smoother surface versus the competition.
- A smooth surface results in a less porous, longer lasting restoration

“Plaque has a difficult time sticking to the surface of my restorations because Alpen Polishers create such a smooth surface...”

Shaun Doyle, DDS  
Nuevo, California



### AVAILABLE GRITS

#### Composite

- F = Fine : 45  $\mu\text{m}$
- UF = Ultra Fine : 25  $\mu\text{m}$

#### Composite PLUS

- F = Fine : 45  $\mu\text{m}$
- UF = Ultra Fine : 15  $\mu\text{m}$
- M = Medium : 110  $\mu\text{m}$

#### Ceramic

- M = Medium : 110  $\mu\text{m}$
- F = Fine : 45  $\mu\text{m}$

#### Ceramic PLUS

- M = Medium : 110  $\mu\text{m}$
- F = Fine : 45  $\mu\text{m}$
- UF = Ultra Fine : 15  $\mu\text{m}$

# Comprehensive Shapes & Sizes

## COMPOSITE

	MINIPOINT	POINT	CUP	SMALL LENS	MINIPOINT	POINT	CUP	SMALL LENS
Cat No.	RP230055RA	RP230059RA	RP230063RA	RP230069RA	RP230072RA	RP230076RA	RP230080RA	RP230086RA
5pk	FINE				EXTRA FINE			
	STEP 1: SHAPING AND SMOOTHING				STEP 2: HIGH LUSTER POLISHING			

## COMPOSITE PLUS

	MINIPOINT	POINT	CUP	SMALL LENS	MINIPOINT	POINT	CUP	SMALL LENS	MINIPOINT	POINT	CUP	SMALL LENS
Cat No.	RP230264RA	RP230265RA	RP230267RA	RP230271RA	RP230291RA	RP230293RA	RP230295RA	RP230297RA	RP230299RA	RP230301RA	RP230303RA	RP230305RA
5pk	MEDIUM				FINE				ULTRA FINE			
	ONE STEP: MIRROR FINISHING				STEP 1: SHAPING AND SMOOTHING				STEP 2: HIGH LUSTER POLISHING			

## CERAMIC

	POINT	CUP	LENS	WHEEL	POINT	CUP	LENS	WHEEL
Cat No.	RP230114RA	RP230122RA	RP230134UM	RP230136UM	RP230139RA	RP230147RA	RP230159UM	RP230161UM
5pk	MEDIUM				FINE			
	STEP 1: REMOVING AND CONTOURING				STEP 2: HIGH LUSTER POLISHING			

## CERAMIC PLUS

	MINIPOINT	FLAME	CUP	SMALL LENS	WHEEL	DISC	LARGE LENS	MINIPOINT	FLAME	CUP	SMALL LENS
Cat No.	RP230313RA	RP230317RA RP230315HP	RP230319RA	RP230321RA	RP230325RA RP230323HP	RP230327HP	RP230329HP	RP230331RA	RP230335RA RP230333HP	RP230337RA	RP230339RA
5pk	MEDIUM							FINE			
	STEP 1: REMOVING AND SHAPING							STEP 2: SMOOTHING			

	WHEEL	DISC	LARGE LENS	MINIPOINT	FLAME	CUP	SMALL LENS	WHEEL	DISC	LARGE LENS
Cat No.	RP230343RA RP230341HP	RP230345HP	RP230347HP	RP230349RA	RP230353RA RP230351HP	RP230355RA	RP230357RA	RP230361RA RP230359HP	RP230363HP	RP230365HP
5pk	FINE			ULTRA FINE						
	STEP 2: SMOOTHING			STEP 3: HIGH LUSTER POLISHING						

### CAD/CAM Restoration PLUS

- SC = Super Coarse : 150 µm
- M = Medium : 110 µm
- F = Fine : 45 µm

### Brush PLUS

- F = Fine : 45 µm

### Universal

- SC = Super Coarse : 150 µm
- F = Fine : 45 µm

### CAD/CAM RESTORATION PLUS

	LENS	SMALL WHEEL	WHEEL	TAPER	CYLINDER	FLAME	CUP	WHEEL	LARGE LENS	DISC
Cat No.	RP230367HP	RP230369HP	RP230371HP	RP230373HP	RP230375HP	RP230379RA RP230377HP	RP230381RA	RP230383HP	RP230385HP	RP230387HP
5pk	SUPER COARSE STEP 1: REMOVING AND SHAPING					MEDIUM STEP 2: SMOOTHING				

	FLAME	CUP	WHEEL	LARGE LENS	DISC	BRUSH PLUS DIAMOND IMPREGNATED	MANDREL Handpiece
Cat No.	RP230391RA RP230389HP	RP230393RA	RP230395HP	RP230397HP	RP230399HP	CUP RP232004RA FLAME RP232006RA	Right Angle RP255047HP RP255049RA
5pk	FINE STEP 3: HIGH LUSTER POLISHING					FINE POLISHES CERAMIC, COMPOSITES AND METALS	FINE STAINLESS STEEL MANDRELS FOR UNMOUNTED POLISHERS

### UNIVERSAL (COMPOSITE AND CERAMIC)

	WIDE CUP	SMALL LENS	LENS	SMALL WHEEL	WHEEL	MINIPOINT	FLAME	BULLET	CUP	WIDE CUP
Cat No.	RP230180UM	RP230186UM	RP230188UM	RP230194UM	RP230196UM	RP230201RA	RP230203RA RP230205HP	RP230207RA RP230209HP	RP230211RA RP230213HP	RP230215RA
5pk	SUPER COARSE ONE STEP: REMOVING, SHAPING AND POLISHING					FINE ONE STEP: REMOVING, SHAPING AND POLISHING				



Unique Easy Store™ Packaging



[www.AlpenRotary.com](http://www.AlpenRotary.com)

# Overview Dentistry/Laboratory Polishers

	Standard	Premium
Working time	Medium / Fast	Very fast
Surface result	Shine	Mirror finish



Category	Non-Diamond		Diamond - the PLUS Line	
	Shape	Grit Size	Shape	Grit Size
Composite, Compomere	MiniPoint	F, UF	MiniPoint	F, EF, UF
	Point	F, UF	Point	F, EF, UF
	Cup	F, UF	Cup	F, EF, UF
	Small Lens	F, UF	Small Lens	F, EF, UF
Ceramic	Point	M, F	Cup	M, F, UF
	Cup	M, F	Disc	M, F, UF
	Wheel	M, F	Flame	M, F, UF
	Lens	M, F	Lens	UF
			Small Lens	M, F
			Large Lens	M, F, UF
			MiniPoint	M, F, UF
			Wheel	M, F, UF
CAD/CAM Restoration (Oxide Ceramic, Zirconia)			Cup	EF, UF
			Cylinder	SC
			Disc	EF, UF
			Flame	EF, UF
			Lens	SC
			Large Lens	EF, UF
			Taper	SC
			Wheel	SC, EF, UF
			Small Wheel	SC
Universal	Bullet	F		
	Cup	F		
	Very Wide Cup	SC		
	Wide Cup	F		
	Flame	F		
	Lens	SC		
	MiniPoint	F		
	Wheel	SC		
Brush	Cup			
	Flame			

**Coltène/Whaledent Inc.**  
 235 Ascot Parkway  
 Cuyahoga Falls, OH 44223 / USA  
 Tel. USA & Canada + 1 800 221 3046  
 + 1 330 916 8800  
 Fax +1 330 916 7077  
 info@coltenewhaledent.com

**Coltène/Whaledent AG**  
 Feldwiesenstrasse 20  
 9450 Altstätten / Switzerland  
 Tel. +41 (0)71 757 53 00  
 Fax +41 (0)71 757 53 01  
 info@coltenewhaledent.ch

Authorized EU Representative:  
**Coltène/Whaledent GmbH + Co. KG**  
 Raiffeissenstrasse 30  
 89129 Langenau / Germany  
 Tel. +49 (0)7345 805 0  
 Fax +49 (0)7345 805 201  
 info@coltenewhaledent.de



P/N 30410C

## Brasseler\* to Alpen Polisher Conversion Chart

Brasseler	Description	Part Number	Shape	Alpen Polisher
ComposiPro™ Universal	0532.21.050	5009095U0	Bullet	RP230207RA
ComposiPro™ Universal	0535.21.060	5009096U0	Cup	RP230211RA
ComposiPro™ Universal	W1.21	5000755U0	Flame	RP230203RA
ComposiPro™ Universal	W5.21	5000767U0	Cup	RP230215RA
Mandrels	303.21	003497U0	Mandrel	RP255049RA
Mandrels	303.11	003494U0	Mandrel	RP255047HP
CeramiPro™ Universal	0301.0	5000600U0	Wheel	RP230196UM
CeramiPro™ Universal	0310.11	5000602U0	Knife Edge	RP230136UM
CeramiPro™ Universal	0311.00	5000603U0	Knife Edge	RP230161UM
CeramiPro™ Universal	0306.00	5000601U0	Wheel	RP230188UM
CeramiPro™ Universal	0316.00	5000604U0	Knife Edge	RP230134UM
CeramiPro™ Universal	0317.00	5000605U0	Knife Edge	RP230159UM
CeramiPro™ Universal	0321.00	5000606U0	Wheel	RP230194UM
CeramiPro™ Universal	0322.00	5000607U0	Knife Edge	RP230186UM
CeramiPro™ Universal	0330.11	5000608U0	Flame	RP230377HP
CeramiPro™ Universal	0340.11	5000611U0	Flame	RP230389HP
CeramiPro™ Universal	0352C.21	5000806U0	Bullet	RP230313RA
CeramiPro™ Universal	0355C.21	5000808U0	Cup	RP230122RA
CeramiPro™ Universal	0362C.21	5000810U0	Bullet	RP230331RA
CeramiPro™ Universal	0365C.21	5000812U0	Cup	RP230147RA
CeramiPro™ Universal	0372C.21.050	5000814U0	Points	RP230349RA
ComposiPro™ One-Step Polishers	C801001.21.060	5008549U11	Cup	RP230267RA
ComposiPro™ One-Step Polishers	C801002.21.05	5008550U1	Points	RP230265RA
ComposiPro™ One-Step Polishers	C801003.31.100	5008551U1	Knife Edge	RP230271RA
ComposiPro™ One-Step Polishers	C801004.31.030	5009222U0	Points	RP230264RA
ComposiPro™ One-Step Brushes	D209002.21.050	5008547U1	Brush Cup	RP232004RA
ComposiPro™ One-Step Brushes	D209003.21.047	5008548U1	Brush Flame	RP232006RA
ComposiPro™ Diacomp	DC1.21	5000702U0	Points	RP230301RA
ComposiPro™ Diacomp	DC1M.21	5000704U0	Points	RP230293RA
ComposiPro™ Diacomp	DC2.21	5000711U0	Knife Edge	RP230305RA
ComposiPro™ Diacomp	DC2M.21	5000713U0	Knife Edge	RP230297RA
ComposiPro™ Diacomp	DC3.21	5000714U0	Cup	RP230303RA
ComposiPro™ Diacomp	DC3M.21	5000716U0	Cup	RP230295RA
CeramiPro™ Dialite®	H2D11.040	5000733U0	Flame	RP230351HP
CeramiPro™ Dialite®	H2DG11.040	5000731U0	Flame	RP230315HP
CeramiPro™ Dialite®	H2DM11.040	5000732U0	Flame	RP230333HP
CeramiPro™ Dialite®	L26DBC.11.260	5000741U0	Knife Edge	RP230327HP
CeramiPro™ Dialite®	L26DGF.11.260	5000742U0	Knife Edge	RP230363HP
CeramiPro™ Dialite®	L26DRM.11.260	5000743UF	Knife Edge	RP230345HP
CeramiPro™ Dialite®	R17D.11.170	5000564U0	Wheel	RP230359HP
CeramiPro™ Dialite®	R17DG.11.170	5000747U0	Wheel	RP230323HP
CeramiPro™ Dialite®	R17DM.11.170	5000748U0	Wheel	RP230341HP
CeramiPro™ Dialite®	W16D.21	5000762U0	Flame	RP230317RA
CeramiPro™ Dialite®	W16DG.21	5000760U0	Flame	RP230335RA
CeramiPro™ Dialite®	W16DM.21	5000761U0	Flame	RP230353RA
CeramiPro™ Dialite®	W17D.21	5000766U0	Cup	RP230355RA
CeramiPro™ Dialite®	W17DG.21	5000764U0	Cup	RP230319RA
CeramiPro™ Dialite®	W17DM.21	5000765U0	Cup	RP230337RA
CeramiPro™ Dialite®	W18D.21	5006778U0	Knife Edge	RP230357RA
CeramiPro™ Dialite®	W18DG.21	5006776U0	Knife Edge	RP230321RA
CeramiPro™ Dialite®	W18DM.21	5006777U0	Knife Edge	RP230339RA
ET Illustra Composite Polisher	132RA	5012646U0	Minipoint	RP230055RA
ET Illustra Composite Polisher	135RA	5012647U0	Cup	RP230063RA
ET Illustra Composite Polisher	138RA	5012648U0	small lens	RP230069RA
ET Illustra Composite Polisher	232RA	5012649U0	Minipoint	RP230072RA
ET Illustra Composite Polisher	235RA	5012640U0	Cup	RP230080RA
ET Illustra Composite Polisher	238RA	5012651U0	small lens	RP230086RA

\*Brasseler and its products are not registered trademarks of Coltène/Whaledent Inc.